

Product Specifications

CL06.2001

INDEX

Tashiro indicator

Tashiro indicator
Tashiro indicateur
Tashiro Indikator
Tashiro indicatore
Tashiro indicador

Titrimetric indicator (NH₄⁺)
1 g C₁₅H₁₄N₃O₂Na + 1 g C₁₆H₁₈ClN₃S·3H₂O / l C₂H₅OH

Mol weight: -
Density: 0,79 g/ml
HS Nr: 38229000
CAS Nr:
UN: 1993
ADR: 3,II



Specifications: