

Product Specifications

CL00.0606**Phosphoric acid 85% v.p.**

Acido fosforico 85% e.p.
Phosphorsäure 85% Rst.
Acido fosfórico 85% t.p.
Acide phosphorique 85% t.p.
Fosforzuur 85% z.z.

**For laboratory use, DAB, Ph. Eur., BP, NF, Ph. Franç., FCC
85+% H₃PO₄**

Mol weight: 98.00 g/mol
Density: 1,71 g/ml
HS Nr: 28092000
CAS Nr: 7664-38-2
UN: 1805
ADR: 8,III

**Specifications:**

Assay	>85%
Chloride (Cl)	<0.0005%
Fluoride (F)	<0.0005%
Nitrate (NO ₃)	<0.0005%
Phosphorous Acid	<0.003%
Sulfate (SO ₄)	<0.005%
Arsenic (As)	<0.0002%
Copper (Cu)	<0.001%
Iron (Fe)	<0.001%
Lead (Pb)	<0.0005%
Zinc (Zn)	<0.001%