

Product Specifications

CL00.1155

ANALYT

di-Potassium hydrogen phosphate a.r.

di-Kaliumwaterstoffosfaat p.a.
di-Potassium hydrogène phosphate p.a.
di-Kaliumwasserstoffphosphat z.A.
di-Potasio hidrógeno fosfato p.a.
di-Potasio idrogeno fosfato p.a.

Mol weight: 174.18 g/mol
Density: 2.44 g/cm³
HS Nr: 28352400
CAS Nr: 7758-11-4
UN:
ADR:

For laboratory use, Ph. Eur.
99+% K₂HPO₄

Specifications:

Assay	>99%
Loss on drying	<1%
Nitrogen (N)	<0.001%
Sodium (Na)	<0.5%
Chloride (Cl)	<0.003%
Sulfate (SO ₄)	<0.005%
Iron (Fe)	<0.001%
Calcium (Ca)	<0.003%
Magnesium (Mg)	<0.003%
Cadmium (Cd)	<0.0001%
Copper (Cu)	<0.0001%
Chromium (Cr)	<0.0001%
Nickel (Ni)	<0.0001%
Heavy Metals as Lead (Pb)	<0.0005%
Arsenic (As)	<0.0001%
Reducing Substances	: passes test