

Product Specifications

CL00.2614

ANALYT

Silver nitrate a.r.

Argent nitrate p.a.
Zilvernitraat p.a.
Argento nitrato p.a.
Plata nitrato p.a.
Silberniträt z.A.
Prata nitrato p.a.
For laboratory use, ACS, ISO, Ph. Eur.
99.99+% AgNO₃

Mol weight: 169.87 g/mol
Density: 4,35 g/cm³
HS Nr: 28432100
CAS Nr: 7761-88-8
UN: 1493
ADR: 5.1,II

**Specifications:**

Assay	>99.99%
Cadmium (Cd)	<0.0001%
Copper (Cu)	<0.0002%
Iron (Fe)	<0.0002%
Manganese (Mn)	<0.0005%
Lead (Pb)	<0.001%
Zinc (Zn)	<0.0001%
Chloride (Cl)	<0.0001%
Nitrite (NO ₂)	<0.05%
Sulfate (SO ₄)	<0.002%
Clarity of the Solution	: Passes Test
Free Acid	: Passes Test
Substances not Precipitated by Hydrochloric Acid	<0.01%