

Product Specifications

CL00.2637

ANALYT

Sulfuric acid 98% a.r.

Acide sulfurique 98% p.a.
Schwefelsäure 98% z.A.
Acido solforico 98% p.a.
Acido sulfúrico 98% p.a.
Zwavelzuur 98% p.a.

For nitrogen det.
97+% H₂SO₄

Mol weight: 98.08 g/mol
Density: 1,84 g/ml
HS Nr: 28070000
CAS Nr: 7664-93-9
UN: 1830
ADR: 8,II



Specifications:

Assay	>97%
Residue after Ignition	<2ppm
Nitrogen (N)	<1ppm