

CL03.0122

HYDREX

Ammonium chloride-Ammonia buffer pH 10 - ISE buffer

Ammoniumchloride-Ammoniak buffer pH 10
Tampon Ammonium chlorure-Ammoniaque pH 10
Ammoniumchlorid-Ammoniak Puffer pH 10
Tampón Amonio cloruro-Amónico pH 10
Tampone Ammonio cloruro-Ammoniacca pH 10

For laboratory use, Ttrations with ISE Cu electrode
54 g NH₄Cl + 350 ml NH₄OH 25% / l H₂O

Mol weight:
Density: 0,99 g/ml
HS Nr: 38229000
CAS Nr:
UN: 2672
ADR: 8,III



Specifications: