

Product Specifications

CL03.0207

HYDREX

Buffer pH 12.00

Buffer pH 12.00
Puffer pH 12.00
Tampon pH 12.00
Tampón pH 12.00
Tampão pH 12.00
Tampone pH 12.00

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 12.00 ($\pm 0.02/20^{\circ}\text{C}$) Phosphate-NaOH**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:
Density: 1,01 g/ml
HS Nr: 38229000
CAS Nr:
UN:
ADR:

Specifications:

pH	= 12.58 @ 0°C
pH	= 12.26 @ 10°C
pH	= 12.00 @ 20°C
pH	= 11.88 @ 25°C
pH	= 11.75 @ 30°C
pH	= 11.53 @ 40°C
pH	= 11.31 @ 50°C
pH	= 11.09 @ 60°C
pH	= 10.88 @ 70°C
pH	= 10.68 @ 80°C
pH	= 10.48 @ 90°C