

Product Specifications

CL03.0215

HYDREX

Buffer pH 6.00

Tampon pH 6.00
Tampão pH 6.00
Buffer pH 6.00
Tampón pH 6.00
Tampone pH 6.00
Puffer pH 6.00

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 6.00 ($\pm 0.02/20^{\circ}\text{C}$) Citrate-NaOH**

**Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:
Density: 1,01 g/ml
HS Nr: 38229000
CAS Nr:
UN:
ADR:

Specifications:

pH	= 6.03 @ 0°C
pH	= 6.01 @ 10°C
pH	= 6.00 @ 20°C
pH	= 6.01 @ 25°C
pH	= 6.02 @ 30°C
pH	= 6.04 @ 40°C
pH	= 6.06 @ 50°C
pH	= 6.09 @ 60°C
pH	= 6.13 @ 70°C
pH	= 6.18 @ 80°C
pH	= 6.24 @ 90°C