

Product Specifications

CL03.0218

HYDREX

Buffer pH 9.00

Tampão pH 9.00

Buffer pH 9.00

Tampone pH 9.00

Tampón pH 9.00

Puffer pH 9.00

Tampon pH 9.00

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 9.00 ($\pm 0.02/20^{\circ}\text{C}$) $\text{H}_3\text{BO}_3\text{-KCl-NaOH}$** **Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:

Density: 1,01 g/ml

HS Nr: 38229000

CAS Nr:

UN:

ADR:

Specifications:

pH	= 9.24 @ 0°C
pH	= 9.11 @ 10°C
pH	= 9.00 @ 20°C
pH	= 8.97 @ 25°C
pH	= 8.93 @ 30°C
pH	= 8.86 @ 40°C
pH	= 8.80 @ 50°C
pH	= 8.75 @ 60°C
pH	= 8.71 @ 70°C
pH	= 8.67 @ 80°C
pH	= 8.64 @ 90°C