

Product Specifications

CL03.0711

HYDREX

Buffer pH 9.00 (colored)

Puffer pH 9.00 (farbcodiert)

Tampão pH 9.00 (colorido)

Buffer pH 9.00 (gekleurd)

Tampon pH 9.00 (coloré)

Tampón pH 9.00 (de color)

Tampone pH 9.00 (colorata)

**For laboratory use, pH measurements, traceable to SRM from NIST
pH 9.00 ($\pm 0.02/20^{\circ}\text{C}$) $\text{H}_3\text{BO}_3\text{-KCl-NaOH}$ (green)****Traceability: Reference standards acc. NIST SRM 189c, 188, 185i,
186lg, 186llg, 187e**

Mol weight:

Density: 1,00 g/ml

HS Nr: 38229000

CAS Nr:

UN:

ADR:

Specifications:

pH	= 9.24 @ 0°C
pH	= 9.11 @ 10°C
pH	= 9.00 @ 20°C
pH	= 8.97 @ 25°C
pH	= 8.93 @ 30°C
pH	= 8.86 @ 40°C
pH	= 8.80 @ 50°C
pH	= 8.75 @ 60°C
pH	= 8.71 @ 70°C
pH	= 8.67 @ 80°C
pH	= 8.64 @ 90°C