

Product Specifications

CL05.2652

TITREX

Sulfuric acid 0.02 mol/l

Acido sulfúrico 0.02 mol/l
Schwefelsäure 0.02 mol/l
Acido solforico 0.02 mol/l
Zwavelzuur 0.02 mol/l
Acide sulfurique 0.02 mol/l

For laboratory use, Titrations & TMAR TA-TAC Det.
1.962 g H₂SO₄ / l H₂O = 0.04 N (±0.00008/20°C)
Traceability: Reference standards acc. NIST SRM 84I

Mol weight: 98.08 g/mol
Density: 1,00 g/ml
HS Nr: 28070000
CAS Nr: 7664-93-9
UN:
ADR:

Specifications: