

Product Specifications

CL12.1919

VINEX

Saccharose 30 weight % solution

Sacarosa 30 p/p % solución
Saccharose 30 gew. % oplossing
Saccarosio 30 p/p % soluzione
Saccharose 30 w/w % Lösung
Saccharose 30 poids % solution

For analysis of wines and vinegars - For Brix refractometer calibration.

300 g C₁₂H₂₂O₁₁ / kg H₂O (= 30°Bx at 20°C)
Traceability: Reference standards acc. NIST SRM 17f

Mol weight: 342.30 g/mol
Density: 1,1270 g/ml
HS Nr: 17019990
CAS Nr: 57-50-1
UN:
ADR:

Specifications: