

Product Specifications

CL12.1923

VINEX

Saccharose 55 weight % solution

Saccharose 55 gew. % oplossing
Saccharose 55 poids % solution
Sacarosa 55 p/p % solución
Saccarosio 55 p/p % soluzione
Saccharose 55 w/w % Lösung

For analysis of wines and vinegars - For Brix refractometer calibration.

550 g C₁₂H₂₂O₁₁ / kg H₂O (= 55°Bx at 20°C)
Traceability: Reference standards acc. NIST SRM 17f

Mol weight: 342.30 g/mol
Density: 1,2580 g/ml
HS Nr: 17019990
CAS Nr: 57-50-1
UN:
ADR:

Specifications: